

Shimadzu unveils tool to evaluate thermal stability of nucleic acid pharmaceuticals

April 18, 2023 | News

For improving the efficiency and reliability of analysis data in the pharmaceutical industry



Japanese firm Shimadzu has released the T_m analysis system both in Japan and internationally. This analyses the thermal stability of nucleic acid pharmaceuticals in development using a UV-VIS spectrophotometer. This is the first system to automate the processes from measurement of the T_m value, a standard index of the thermal stability of nucleic acids, to data analysis.

The T_m value is indispensable for research and development in the nucleic acid pharmaceutical market, which is expected to have an annual growth rate of approximately 17% by 2030.

The T_m analysis system automatically performs everything from sample measurements to data analysis in one step. The system consists of a UV-VIS spectrophotometer in the UV-Vis series, the TMSPC-8i 8-series thermoelectrically cooled cell holder, and the LabSolutions UV-Vis T_m measurement analysis software. The newly developed micro cell limits evaporation of the sample, enabling trace quantity measurements, so precious nucleic acid samples are not wasted.

The newly designed software is compliant with standards related to data integrity (data integrity concomitant on measurement) as required in the pharmaceutical industry.